
# EnzyFluo<sup>TM</sup> NF<sub>K</sub>B Phosphorylation Assay Kit (ENF<sub>K</sub>B-100)

Fluorimetric Cell-Based Assay for p65/RelA (S536) Phosphorylation Status

## DESCRIPTION

Nuclear factor-kappa B (NF $\kappa$ B) is a transcription factor that plays a central role in many physiological processes, e.g. inflammation, tumorigenesis, and apoptosis. NF $\kappa$ B is activated by a wide variety of stimuli, including inflammatory cytokines such as TNF- $\alpha$ . NF $\kappa$ B is a dimer composed of members of the Rel family of proteins: p65/RelA, c-Rel, RelB, NF $\kappa$ B1/p50, and NF $\kappa$ B2/p52. Phosphorylation of p65/RelA at Ser-536 results in decreased nuclear export and enhanced p65/RelA-dependent transcription. BioAssay Systems' cell-based ELISA measures phosphorylated p65(S536) (pNF $\kappa$ B) in whole cells and normalizes the signal to the total protein content. This simple and efficient assay eliminates the need for cell lysate preparation and can be used to study NF $\kappa$ B regulation in short-term and long-term assays.



## **KEY FEATURES**

*New and Improved.* Total assay time reduced from typical 9 hours to 5 hours.

Simple. No cell lysis necessary. Cells can be directly cultured in 96-well plates.

Convenient. Total protein and  $pNF\kappa B$  can be measured in the same sample.

## APPLICATIONS

Determination of NF $\kappa$ B p65 (S536) phosphorylation status in whole cells. *Evaluation* of direct and indirect modulation of NF $\kappa$ B p65 phosphorylation. *Species tested*: human, mouse.

## **KIT CONTENTS**

10× Wash Buffer:	25 mL	Blocking Buffer:	25 mL
Protein Stain:	6 mL	HRP Substrate:	6 mL
pNFκB-Ab1:	10 µL	HRP-Ab2	10 µL

**Storage conditions:** This kit is shipped on ice. Upon delivery, store all reagents at -20°C. Shelf life of 6 months after receipt.

**Precautions:** reagents are for research use only. Normal precautions for laboratory reagents should be exercised while using the reagents. Please refer to material safety data sheet for detailed information.

### ASSAY PROCEDURE

Important:

- To avoid cross-contamination, change pipette tips between additions of each reagent or sample. We recommend the use of a multi-channel pipette. Use separate reservoirs for each reagent. Prior to Assay, dilute 10× Wash Buffer in dH<sub>2</sub>O to prepare 250 mL 1× Wash Buffer.
- 2. It is recommended that samples be assayed in triplicate or higher.
- 3. Two different blanks are necessary. For each *plate* include a Protein Blank (no cells) in triplicate. For each *sample* include a Sample Blank (cells w/ only Ab2) in triplicate. The blanks are used to determine background fluorescence for total protein and pNFκB respectively.

#### A. Culture and Treat Cells

1. Seed 100  $\mu$ L of 1-3 × 10<sup>4</sup> adherent cells (or 4-10 × 10<sup>4</sup> suspension cells) into each well of a black 96-well culture plate. Add 100  $\mu$ L of culture media without cells into three wells for the Protein Blank. Incubate overnight at 37°C in a cell culture incubator.

Note: The cell number to be used depends on the cell line and NF  $\kappa\!B$  phosphorylation status.

- 2. Treat the cells as desired (e.g. with ligands or drugs).
- 3. Prepare formaldehyde solutions (*warning*: formaldehyde is toxic. Use chemical hood and wear appropriate gloves and eye protection):

For adherent cells, prepare 4% formaldehyde by mixing 1.3 mL of 37% formaldehyde and 10.7 mL of 1× Wash buffer. Simply fix cells in each well by replacing the medium with 100  $\mu$ L of 4% formaldehyde.

For suspension cells, prepare 8% formaldehyde by mixing 2.6 mL of 37% formaldehyde and 9.4 mL of 1× Wash buffer. Centrifuge the plate at 500g for 15 min at 4°C and carefully remove as much media as possible without disturbing the cell pellet (repeat this for suspension cells with each wash step below). Fix the cells in each well by adding 100  $\mu$ L of 8% formaldehyde to cell pellet.

Cover the plate and incubate for 20 min at room temperature. Alternatively, the plate containing the fixed cells can be sealed and stored for up to 2 weeks at  $2-8^{\circ}$ C.

- 4. Remove the formaldehyde solution and wash the cells 3 times with 150  $\mu L$  of 1× Wash Buffer. Each wash step should be performed for 1 min with gentle shaking.
- 5. Prepare Quench Buffer by mixing 2.2 mL of 3%  $H_2O_2$  and 8.8 mL of 1× Wash Buffer.

Remove the Wash Buffer and add 100  $\mu L$  of Quench Buffer to each assay well. Cover plate and incubate for 20 min at room temperature.

- 6. Remove the Quench Buffer and wash the cells 3 times with 150  $\mu L$  of 1× Wash Buffer.
- 7. Remove the Wash Buffer, and add 100  $\mu$ L of Blocking Buffer. Cover plate and incubate for 1 hr at room temperature.

#### B. Add Primary Antibodies (Ab1)

- 1. Prepare 55  $\mu$ L of primary antibody Ab1 Mixture for each well by mixing pNF $\kappa$ B -Ab1 into Blocking Buffer in a 1:625 dilution.
- Remove the Wash Buffer from all assay wells. Add 50 µL of the Blocking Buffer to the Sample Blank wells and 50 µL of Ab1 Mixture to the Sample wells. Cover plate and incubate for 90 min at room temperature or overnight at 2-8°C with gentle shaking.
- Remove the Ab1 Mixture and wash the cells 3 times with 150 µL of 1× Wash Buffer. Each wash step should be performed for 1 min with gentle shaking.

#### D. Add Secondary Antibodies (Ab2)

- 1. Prepare 55  $\mu$ L of secondary antibody Ab2 Mixture for each well by mixing HRP-Ab2 into Blocking Buffer in a 1:625 dilution.
- 2. Remove Wash Buffer and add 50  $\mu L$  of the Ab2 Mixture to all assay wells. Cover plate and incubate for 90 min at room temperature with gentle shaking.

#### E. Detection

- Remove the Ab2 Mixture from each well and thoroughly wash the cells 5 times with 150 µL of 1× Wash Buffer. Each wash step should be performed for 1 min with gentle shaking.
- 2. Immediately before use, add 6  $\mu L$  3%  $H_2O_2$  to the provided 6 mL HRP Substrate (for partial plate assay, adjust the volumes accordingly). Remove the Wash Buffer from the plate and add 50  $\mu L$  of mixed HRP Substrate to each well. Incubate for 30 min at room temperature in the dark.
- 3. Add 50 µL of Protein Stain to each well and incubate for an additional 5 min at room temperature in the dark.
- 4. Read the plate at  $\lambda_{\text{ex/em}}$  = 530/585 nm for phosphorylated NF $\kappa$ B (pNF $\kappa$ B) and at  $\lambda_{\text{ex/em}}$  =360/450 nm for total protein.

## CALCULATION

Calculate the mean pNF<sub>K</sub>B fluorescence intensities at 530/585nm for the Sample Blank wells ( $F_{BLK}$  <sub>pNF<sub>K</sub>B</sub>) and Sample wells ( $F_{SAMPLE}$  <sub>pNF<sub>K</sub>B</sub>). Also calculate the mean protein fluorescence intensities at 360/450nm for the Protein Blank (no cells) well ( $F_{BLK}$  <sub>prot</sub>) and Sample wells ( $F_{SAMPLE}$  <sub>Prot</sub>). Subtract the mean pNF<sub>K</sub>B fluorescence of the Sample Blank wells from the pNF<sub>K</sub>B fluorescence of the Sample wells to yield  $\Delta F$  values for the pNF<sub>K</sub>B ( $\Delta F_{pNF_KB}$ ). Subtract the mean protein fluorescence of the Protein Blank (no cells) well from the protein fluorescence of the Sample wells to yield  $\Delta F$  values for the pNF<sub>K</sub>B fluorescence value of the Sample wells to yield  $\Delta F$  values for the total protein ( $\Delta F_{PROT}$ ).

$$\Delta \overline{F}_{pNFkB} = \overline{F}_{pNFkB}^{SAMPLE} - \overline{F}_{pNFkB}^{BLK}; \quad \Delta \overline{F}_{Prot} = \overline{F}_{Prot}^{SAMPLE} - \overline{F}_{Prot}^{BLK}$$


Normalized phosphorylated NF<sub>K</sub>B is calculated as,

Normalized pNF<sub>K</sub>B = 
$$\frac{\Delta F_{PNFKB} / \Delta F_{Prot}}{(\Delta \overline{F}_{PNFKB} / \Delta \overline{F}_{Prot})_{o}}$$

where  $(\Delta F_{pNF_KB} / \Delta F_{Prot})_o$  is the control reference value (e.g. time zero in kinetic studies or untreated wells in drug potency studies.)

## MATERIALS REQUIRED BUT NOT PROVIDED


37% formaldehyde (Sigma, cat # F8775); 3% H<sub>2</sub>O<sub>2</sub> (Sigma, cat # 323381); black cell culture 96-well plate: available separately at BioAssay System (cat# P96BCC) or at VWR (cat# 82050-748); plate sealers: available separately at BioAssay Systems (cat# AB96SL) or at Sigma (cat# A5596); deionized or distilled water; pipetting devices; cell culture incubators; centrifuge tubes; fluorescence plate reader capable of reading at  $\lambda_{ex/em} =$ 530/585 nm and at  $\lambda_{ex/em} =$  360/450 nm.



*Left*: Phosphorylation of NF $\kappa$ B p65(S536) in human breast cancer cell line MDA-MB-231 and murine NIH-3T3 fibroblast cells after stimulation with human TNF- $\alpha$ . *Right*: Kinetics of NF $\kappa$ B p65(S536) phosphorylation in MDA-MB-231 cells after TNF- $\alpha$  stimulation.

## LITERATURE

- Neumann M and Naumann M (2007). Beyond IkappaBs: alternative regulation of NF-kappaB activity. FASEB J. 21(11):2642-54
- Jiang X et al (2003). The NF-kappa B activation in lymphotoxin beta receptor signaling depends on the phosphorylation of p65 at serine 536. J Biol Chem. 278(2):919-26.
- Gutierrez H et al (2008). Nuclear factor kappa B signaling either stimulates or inhibits neurite growth depending on the phosphorylation status of p65/ReIA. J Neurosci. 28(33):8246-56.



# EXAMPLE OF A 96-WELL ASSAY PLATE LAY-OUT