

Safety Data Sheet

SECTION 1: Identification of the substance/mixture and of the company/undertaking

Identification of the substance or mixture

Product code 25200072
Product name Trypsin, 0.25% (1X) with EDTA 4Na, liquid

Company/undertaking identification

Life Technologies Corporation
5781 Van Allen Way
PO Box 6482
Carlsbad, CA 92008
+1 760 603 7200

Life Technologies
5250 Mainway Drive
Burlington, ONT
CANADA L7L 6A4
800/263-6236

24 hour Emergency Response for Hazardous Materials [or Dangerous Goods] Incident. Spill, Leak, Fire, Exposure, or Accident. Call CHEMTREC Within the USA + Canada: 1-800-424-9300 and +1 703-527-3887
Outside the USA + Canada: +1 703-741-5970

Country specific Emergency Number (if available):

CHEMTREC Brazil (Rio De Janeiro) +(55)-2139581449 (português)

CAUTION: For use as a raw material component in further manufacturing applications.

SECTION 2: Hazards identification

GHS - Classification

Signal Word
DANGER

Hazard pictograms



Health hazards

Respiratory sensitization

Category 1

Physical hazards

Not Hazardous

Revision date 23-Aug-2017
Product code 25200072

Page 1 / 7
Product name Trypsin, 0.25% (1X) with EDTA 4Na, liquid

Environmental hazards

Not Hazardous

Hazard Statements

H334 - May cause allergy or asthma symptoms or breathing difficulties if inhaled

Precautionary Statements**Prevention**

P261 - Avoid breathing dust/fume/gas/mist/vapors/spray

P284 - In case of inadequate ventilation wear respiratory protection

Response

P304 + P340 - IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing

P342 + P311 - If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician

Storage

Not Applicable

Disposal

P501 - Dispose of contents/ container to an approved waste disposal plant

Other hazards

Not Applicable

HMIS

Health	4
Flammability	0
Reactivity	0

SECTION 3: Composition/information on ingredients

Component	CAS-No	EINECS-No	Weight %
Trypsin 9002-07-7 (0.1-1.0)	9002-07-7	232-650-8	0.1-1.0

We recommend handling all chemicals with caution.

SECTION 4: First aid measures

Description of first aid measures

Skin contact	Rinse with plenty of water . Immediate medical attention is not required.
Eye contact	Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do
Ingestion	Not expected to present a significant ingestion hazard under anticipated conditions of normal use. If you feel unwell, seek medical advice.
Inhalation	Not expected to be an inhalation hazard under anticipated conditions of normal use of this material. Consult a physician if necessary.
Notes to Physician	Treat symptomatically.

Most important symptoms and effects, both acute and delayed

H334 - May cause allergy or asthma symptoms or breathing difficulties if inhaled

Indication of any immediate medical attention and special treatment needed

IF INHALED: Remove person to fresh air and keep comfortable for breathing.

SECTION 5: Firefighting measures

Extinguishing media

Suitable extinguishing media	Water spray. Carbon dioxide (CO ₂). Foam. Dry chemical.
Unsuitable extinguishing media	No information available.

Special hazards arising from the substance or mixture Not known.

Advice for fire-fighters Standard procedure for chemical fires.

SECTION 6: Accidental release measures

Personal precautions, protective equipment and emergency procedures

Ensure adequate ventilation. Always wear recommended Personal Protective Equipment. Use personal protection equipment. See Section 8 for more detail.

Environmental precautions

No special environmental precautions required. Avoid discharge into drains and waterways whenever possible.

Methods and material for containment and cleaning up

Soak up with inert absorbent material.

Reference to other sections

See section 8 for more information.

SECTION 7: Handling and storage

Precautions for safe handling

Use personal protective equipment as required. No special handling advices are necessary.

Conditions for safe storage, including any incompatibilities

Keep in a dry, cool and well-ventilated place. Keep in properly labeled containers.

Specific end use(s)

CAUTION: For use as a raw material component in further manufacturing applications.

Revision date 23-Aug-2017
Product code 25200072

Page 3 / 7
Product name Trypsin, 0.25% (1X) with EDTA 4Na, liquid

SECTION 8: Exposure controls/personal protection

Control parameters

Chemical Name	OSHA PEL	OSHA PEL (Ceiling)	ACGIH OEL (TWA)	ACGIH OEL (STEL)
Trypsin	None	None	None	None

Engineering measures Ensure adequate ventilation, especially in confined areas

Exposure controls

Personal Protective Equipment

Respiratory protection In case of insufficient ventilation wear respirators and components tested and approved under appropriate government standards.

Hand protection Wear suitable gloves Glove material: Compatible chemical-resistant gloves.

Eye protection Tight sealing safety goggles

Skin and Body Protection Wear suitable protective clothing

Hygiene measures Handle in accordance with good industrial hygiene and safety practice

Environmental exposure controls

No special environmental precautions required.

SECTION 9: Physical and chemical properties

Information on basic physical and chemical properties

Appearance	liquid	
Odor	No data available	
Melting point / melting range	°C Mixture has not been tested	°F Mixture has not been tested
Boiling point / boiling range	°C Mixture has not been tested	°F Mixture has not been tested
Flash point	°C Mixture has not been tested	°F Mixture has not been tested
Autoignition Temperature	°C Mixture has not been tested	°F Mixture has not been tested
Decomposition temperature	°C Mixture has not been tested	°F Mixture has not been tested
Evaporation rate	No data available	
Flammability (solid, gas)	No data available	
Upper explosion limit	Mixture has not been tested	
Lower explosion limit	Mixture has not been tested	
Vapor Pressure	Mixture has not been tested	
Relative density	Mixture has not been tested	
Specific gravity	No data available	
Solubility	No data available	
Partition coefficient: n-octanol/water	No data available	
Explosive properties	Mixture has not been tested	
Other information	No data available	

SECTION 10: Stability and reactivity

Reactivity	None known.
Chemical stability	Stable under normal conditions.
Possibility of hazardous reactions	Hazardous reaction has not been reported.
Conditions to avoid	No information available.
Incompatible materials	No dangerous reaction known under conditions of normal use.
Hazardous decomposition products	No data available.

SECTION 11: Toxicological information

Information on toxicological effects

Chemical Name	LD50 (oral, rat/mouse)	LD50 (dermal, rat/rabbit)	LC50 (inhalation, rat/mouse)
Trypsin	> 5 g/kg (Rat)	No data available	No data available

Principal Routes of Exposure

Irritation	Conclusive but not sufficient for classification
Corrosivity	Conclusive but not sufficient for classification
Sensitization	May cause sensitization by inhalation
STOT - Single Exposure	Conclusive but not sufficient for classification
STOT - Repeated Exposure	Conclusive but not sufficient for classification
Carcinogenicity	Conclusive but not sufficient for classification
Mutagenicity	Conclusive but not sufficient for classification
Reproductive toxicity	Conclusive but not sufficient for classification
Aspiration hazard	Conclusive but not sufficient for classification

SECTION 12: Ecological information

Toxicity

The environmental impact of this product has not been fully investigated.

Chemical Name	Freshwater Algae Data	Water Flea Data	Freshwater Fish Species Data	Microtox Data	log Pow
Trypsin	No data available	No data available	No data available	No data available	No data available

Persistence and degradability No information available.

Bioaccumulative potential No information available.

Results of PBT and vPvB assessment

No information available.

Other adverse effects No information available.

SECTION 13: Disposal considerations

Waste treatment methods

The generation of waste should be avoided or minimized wherever possible. Empty containers or liners may retain some product residues. This material and its container must be disposed of in accordance with approved disposal technique. Disposal of this product, its solutions or of any by-products, shall comply with the requirements of all applicable local, regional or national/federal regulations.

SECTION 14: Transport information

IATA / ADR / DOT-US / IMDG

Not regulated in the meaning of transport regulations.

UN number Not Applicable

UN proper shipping name Not Applicable

Transport hazard class(es) Not Applicable

Packing group Not Applicable

Environmental hazards Not Applicable

Special precautions for user Not Applicable

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not Applicable.

SECTION 15: Regulatory information

Component	US TSCA
Trypsin 9002-07-7 (0.1-1.0)	Listed

US Federal Regulations

SARA 313

This product is not regulated by SARA.

Clean Air Act, Section 112 Hazardous Air Pollutants (HAPs) (see 40 CFR 61)

This product does not contain HAPs.

US State Regulations

California Proposition 65

This product does not contain any Proposition 65 chemicals.

WHMIS Hazard Class

D2A - Very toxic materials



This product has been classified in accordance with the hazard criteria of the Controlled Products Regulations (CPR) and the MSDS contains all the information required by the CPR.

SECTION 16: Other information

Reason for revision SDS sections updated.
Revision number 4
Revision date 23-Aug-2017

CAUTION: For use as a raw material component in further manufacturing applications.

References

- ECHA: <http://echa.europa.eu/>
- TOXNET: <http://toxnet.nlm.nih.gov/>
- eChemPortal: <http://www.echemportal.org/>
- LOLI database: <https://www.chemadvisor.com/loli-database>

"The above information was acquired by diligent search and/or investigation and the recommendations are based on prudent application of professional judgment. The information shall not be taken as being all inclusive and is to be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Since the Company cannot control the actual methods, volumes, or conditions of use, the Company shall not be held liable for any damages or losses resulting from the handling or from contact with the product as described herein. THE INFORMATION IN THIS SDS DOES NOT CONSTITUTE A WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE"

End of Safety Data Sheet

Revision date 23-Aug-2017
Product code 25200072

Page 7 / 7
Product name Trypsin, 0.25% (1X) with EDTA 4Na, liquid